

## NORHOF LN2 Microdosing System #608

### FTIR detector filling system

- 💧 Universal automatic filling system suitable for basically all FTIR detectors
- 💧 We produce adapters tailor-made to work smoothly on systems of Bruker, Netzsch, Shimadzu, Thermo Fisher, Perkin Elmer, Agilent, Horiba, Mettler Toledo and many more
- 💧 Plug and play refilling system, to refill the detector at a pre-set time interval
- 💧 Usable for overnight measurements
- 💧 Possible to control extra filling:
  - ⇒ to prevent an automatic fill action will disturb your automatic measurements
- 💧 Fully safe to use inside the laboratory
- 💧 Always ready to use
- 💧 Very slow and gentle filling rate
- 💧 Easy setup, on most FTIR detectors, no tooling needed



*Norhof auto-fill adapter on a FTIR Dewar*

### Norhof #608 LN2 microdosing systems



Norhof manufactures LN2 microdosing systems. Liquid Nitrogen (LN2) is used as the cooling medium and is taken from a storage vessel (Dewar) with low pressure (max. 300 mBar) and delivered (pumped) through a fill line to the application in a micro dosing way.

The Norhof LN2 microdosing system is designed to overcome the drawbacks of LN2 under pressure in which a solenoid valve is used to switch the supply ON / OFF. You may compare the Norhof system with a water tap, but instead of giving water, it gently gives liquid nitrogen, with an adjustable flow, possible to regulate from some drops, up to 0,5 Liter/minute. Our pump can pump LN2 up to 5 meters above the pump itself

*Norhof #608 pump, mounted on a 35 Liter Dewar*


### Working principle

The pressure above the liquid level inside the Dewar is built by heating a small amount of liquid in the bottom of the Dewar.

With only up to 100 mBar of overpressure, the liquid will gently rise out of the rise pipe and fall into the fill hose.

Because we evaporate some LN2 to build pressure, there is no adding of ice inside the Dewar, such as with manual systems which use air from the environment.

When LN2 is required, a small overpressure is generated by a small heater element in the LN2, and liquid flows out of the system like water from a tap, without spilling, noise, vibrations etc.



## #608 LN2 microdosing

The reservoir Dewar can safely stand next to your working place, ready for use. The #608 pump is adjusted with a certain time interval, which is shorter than the hold-time of your detector. After this interval time is past, the detector is carefully and slowly refilled. Refilling is done with a very gentle steady flow. As soon as the detector is full, the pump stops, and will start again after the set time interval, or when another signal is sent to the pump. In this way you are sure that your detector is always cold.

## #608 Technical Specifications

Flow rate	Up to 0,5 liter per minute		
Average working pressure	100 mBar (adjustable with the supplied software)		
Maximum working pressure	300 mBar		
Reaction time	+/- 1 minutes for cooling down the fill line (with 1.6 meters fill line)		
Power connection	115V / 230V AC with supplied power supply or 12-24 Volt AC/DC		
Power consumption	average 10 watt, during pumping 50 watts		
Model	<b>608-100-035</b>	<b>608-100-050</b>	<b>608-100-100</b>
Storage container volume	35 Liter	50 Liter	100 Liter
Outside dimensions (diameter)	460 mm	460 mm	500 mm
Height dimensions	802 mm	917 mm	1235 mm
Weight (empty, full)	15 / 42 kg	17 / 59 kg	34 / 115 kg
Static holding time (days)	112 days	122 days	164 days
Standard fill line	1,6 m line length, 6,25 mm OD, 4 mm ID, Teflon FEB tube, with 32 mm OD Armaflex foam insulation		
System includes	Dewar, pump, fill line 1,6 m, phase separator, fill tubes, power supply, cables, 1 level sensor, PC software.		
Working modes	Automatic fill control with timer (1 sensor)		
External control	5 volt signals for ON/OFF and RS232 signals for ON/OFF		
PC software	Monitor software, to monitor pump and data logging		
Alarms/warning acoustical/ visual / mechanical	Dewar empty, Dewar 4 liter LN2 left, broken sensor(s), frozen alarm, mechanical overpressure protection valve.		
Options	Tailor-made Adaptor with sensor for all brands of detectors Custom built adaptor to fixate sensor(s) on application Transport trolley 5 wheels (10 cm height) Stand for pump (when Dewar is refilled)		
Dewar Trolley option	900-400-000	900-400-000	900-400-000
Floor Stand option	900-400-100	900-400-100	900-400-100
Custom made adapter	600-400-000	600-400-000	600-400-000

## #608 recommended options:



**Norhof Dewar Trolley**  
for safe and easy transport of the LN2 Dewar



**Norhof Floor Stand**  
for storing the pump in vertical position during refill of the LN2 Dewar

## #608 Series FTIR Adaptors



Standard we have adaptors to fit on systems of Bruker, Netzsch, Shimadzu, Thermo Fisher, Perkin Elmer, Agilent, Horiba, Mettler Toledo but also for other detectors we have adaptors available or make a tailor-made design.

The majority of these system finds use in the field of laboratories, research centers where FTIR measurements are made and where a safe and unattended refill of LN2 is required.

## #608 series applications:


- Unattended automatic refill of LN2 for FTIR detectors;
- Refilling of laser detectors
- Overnight measurements
- Automatic sequences of measurements
- Auto filling of other small LN2 reservoirs, with very low LN2 flowrate



## Software

To display sensor temperatures, vessel pressure, status of LED's on the pump etc. our Norhof Monitoring software is included with any #608 series pump. This software works under Windows '98 - 2000 - ME - NT - Vista - W7 – Windows 10 – Windows 11

With each #608 system a software Datalogger function is included. This recorder allows you to document any parameter value over time and to preset the time interval for your detector.



## #608 advantages:

- **the system is extremely safe**  
the operator is not in direct contact with LN2
- **the system is time saving**  
the operator does not need to pour LN2 several times in a run
- **the system can cool the detector just with a press on the start button**  
this means that the detector is cold 24 hours a day
- **there is no LN2 valve required**  
that implies no unnecessary heat input
- **there is no additional control unit required**  
which adds to a clean and elegant setup
- **there is a very low thermal connection to the ambient temperature**  
this means the system is extremely economic in stand-by modes (typical usage less than 0,5 Liter / day)
- **the system can deliver LN2 liquid with a flow optimized for the application**  
without noise, vibration, excessive waste, etc.
- **the system is prepared to be connected to a PC**  
perfect for monitoring and data logging
- **P.E.D. 99/36/EC (Pressure European Directive) for pressurized vessels does not apply for this system**  
The maximum possible pressure is 300mBar. Therefore this system is allowed to be used inside the lab, near your working place, without danger.